

# Infusomat® P

Broader application spectrum through extended performance features



Compatible with standard i.v. administration sets


Volumetric infusion pump	<b>Classification:</b> CF, defibrillation-proof, protection class I																		
Delivery performance	<b>Delivery range and rates:</b> 0,1 – 999.9 ml/h; smallest increment that can be entered is 0.1 ml/h <b>Accuracy of delivery rate setting:</b> typical $\pm 5\%$ with Intrafix® P i.v. administration set, measurement time 24 h, temperature 22°C																		
Function spectrum	<b>Volume selection:</b> 0.1...9999.9 ml / <b>Time selection:</b> max. 99 h 59 min <b>Automatic calculation of delivery rate:</b> from volume and time settings <b>Bolus administration, adjustable parameters:</b> bolus volumes 0,1...99.9/500.0 ml, bolus rates 1...999.9 ml/h, interval times max 99 h 59 min <b>Dose calculation:</b> automatic calculation of delivery rate, ml/h from input $\mu\text{g}/\text{mg}/\text{IE}$ or mmol per ml volumes, time-based per min/h/day (24 h) and/or body weight- and time-based, per kg and/or per min/h/day (24 h) <b>History function:</b> time-based storage of 350 events for therapy analysis <b>Piggyback function:</b> sequential delivery of primary and secondary infusion by switch-over <b>Standby infusion pause:</b> max. 24 h, set values remain stored <b>Drug selection:</b> up to 9 drug names displayed <b>Deactivatable drop control:</b> linked to volume setting <b>Data lock:</b> keyboard lock prevents alteration of preset values <b>Battery capacity display:</b> in h and min																		
Configuration	<b>Integrated service program for application-specific configuration and technical diagnosis of equipment</b>																		
Information	<b>Large, illuminated display:</b> shows delivery rate, mains/battery operation, "pump in operation", alarm situation, functions/operator guidance and alarm sources are displayed as text messages <b>Info prompt:</b> current volume, current time, dosage, battery capacity, hours of operation																		
Safety system	<b>Drop control:</b> operation possible with or without drop sensor <b>Flow clamp:</b> prevents free flow when pump door is opened <b>Air detector:</b> minimizes risk of air infusion <table border="1"> <tr> <td><b>Selectable occlusion alarm pressure settings:</b></td> <td>Low</td> <td>High</td> </tr> <tr> <td></td> <td>approx. 0.6 bar</td> <td>approx. 1.2 bar</td> </tr> </table> <table border="1"> <tr> <td><b>Alarm response time for rate</b></td> <td>1 ml/h</td> <td>17 min</td> <td>30 min</td> </tr> <tr> <td></td> <td>25 ml/h</td> <td>45 sec</td> <td>72 sec</td> </tr> <tr> <td></td> <td>100 ml/h</td> <td>12 sec</td> <td>18 sec</td> </tr> </table> <b>Bolus volume in system</b> 0.3 ml                      0.5 ml <b>Visual and acoustic alarms:</b> pressure alarm, drop alarm, air alarm, standby alarm, low battery pre-alarm, battery alarm ... <b>Staff call:</b> max. 24V / 1A/ 24VA (VDE 0834) <b>KOR, Keep-Open-Rate:</b> for delivery rates > 10 ml/h = 3 ml/h and < 10 ml/h = 1 ml/h	<b>Selectable occlusion alarm pressure settings:</b>	Low	High		approx. 0.6 bar	approx. 1.2 bar	<b>Alarm response time for rate</b>	1 ml/h	17 min	30 min		25 ml/h	45 sec	72 sec		100 ml/h	12 sec	18 sec
<b>Selectable occlusion alarm pressure settings:</b>	Low	High																	
	approx. 0.6 bar	approx. 1.2 bar																	
<b>Alarm response time for rate</b>	1 ml/h	17 min	30 min																
	25 ml/h	45 sec	72 sec																
	100 ml/h	12 sec	18 sec																
Main/battery operation, dimensions	<b>Rated voltage:</b> 230V, 50/60 Hz ~ (see ordering information) <b>External extra-low voltage:</b> 12 V --- <b>Battery type:</b> NiCd (rechargeable) 6 x 1.2V; 1.2Ah <b>Battery operating/recharge time:</b> >3.0 h at max. delivery rate / > 16 h <b>Dimensions (W×H×D) / weight:</b> 140×240×200 mm/approx. 3.1 kg																		
Interfaces	<b>Interface port:</b> RS 232 for automatic generation of therapy documentation /therapy management <b>MFC plug-in contact for fm intensive workstations:</b> data transmission, staff call, voltage supply and potential equalization																		
Software, updates, system extension	<b>Software upgrades/updates via interface:</b> extended performance spectrum / updating performance parameters using PC programs																		
Warranty, technical inspections	<b>Warranty:</b> 24 months from delivery date <b>Technical inspections:</b> at two-year intervals																		

## The safety system for less stress

- The operator guidance features provide all the information required for safe set up and operation. The main parameters are always in view.
- A flow clamp prevents "free flow" when the pump door is opened. The air detector minimises the risk of air infusions.
- An adjustable shut-down pressure setting reduces the bolus volume and shortens alarm response time in the event of system occlusion.

# therapeutic goals more easily and safely

## Higher performance, more consistent long-term delivery rates



Infusomat® P – even more consistent delivery rates during long-term infusions. Decreasing resilience of PVC infusion lines is compensated electronically. The curve was calculated on an Intrafix® i.v. administration set at a delivery rate of 100 ml/h at 22°C.

■ The optimized pump control of the Infusomat® P provides more consistent delivery rates during long-term infusions. The declining delivery rate typically observed with PVC infusion lines with increasing infusion time is compensated electronically.

■ Improved pump control also allows delivery rates up to 999.9 ml/h.

## New functions for a broader application spectrum

■ Switching-over between a primary and a secondary infusion in "piggyback" mode is easier. The current infusion program is interrupted for dosage of the second infusion.

■ The Infusomat® P now features a "history" function for clinical course analysis. A data memory stores a record of the last 350 events and settings. The history can be transferred to computer via interface.

## Easy setting, automatic dose calculation

■ Direct rate input or automatic calculation of the delivery rate based on volume and time settings. The rate can be adjusted, if necessary, without stopping the infusion.

■ Bolus application by pre-selectable parameters.

■ Dosage units are administered according to body weight and/or time. Easy parameter input mode, automatic calculation of delivery rate.

## More flexibility with standard sets

■ Compatibility of standard i.v. administration sets like Intrafix® P eliminates the need to change the set when switching from gravity to pump infusion.

■ A selection of further standard sets from other manufacturers can be activated in the service program if desired.



Fast and safe fixation of the set. Flow clamp prevents "free flow" when the pump door is opened.



Can be used, either on stands, rails or as a table-top device. Integrateable into fluid manager system and info intensive workstations.

# Infusomat® P – Ordering information

## Infusion pump and accessories

Volumetric infusion pump		Code No.
Infusomat® P	230 V	0871 2379
Infusomat® P	200/240 V	0871 2212
Infusomat® P	100/120 V	0871 2387
<b>Equipment supplied:</b>		
Drip sensor with 2 adapters for drip chamber large/small		
Universal clamp fm, rotatable		
Suggested accessories		Code No.
MFC-connecting lead for potential equalisation		0870 1628
MFC-Connecting lead for staff call		0871 1682
MFC-connecting lead for ambulances (12V)		0871 1674
MFC-interface lead for RS-232, electrically isolated		0871 1661
Universal adapter for drip sensor		0871 1747
Infusion tripod, short with drip chamber holder		0870 1644

## Disposables

Intrafix® P i.v. administration sets	Code No.
Intrafix® P, silicone-coated spike, length of tubing: approx. 180 cm	0407 0038
Intrafix® Air P, length of tubing: approx. 150 cm	0406 2957
Intrafix® Air P, silicone-coated spike, length of tubing: approx. 150 cm	0406 2930
Intrafix® Air P, tubing length approx. 180 cm	0406 2981
Intrafix® Air P, silicone-coated spike, length of tubing: approx. 180 cm	0406 2990
Intrafix® Air RV P, with back-check valve for parallel infusions, length of tubing: approx. 150 cm	0406 3252
Intrafix® Air RV P, with back-check valve for parallel infusions, length of tubing: approx. 180 cm	0406 3287
Intrafix® Air-matic® P, length of tubing: approx. 150 cm	0406 3953
Intrafix® Air-matic® P, length of tubing: approx. 180 cm	0406 3988
Intrafix® Air-matic® P, silicone-coated spike, length of tubing: approx. 180 cm	0406 2973
Intrafix® Air-matic® RV P, with back-check valve for parallel infusions, silicone-coated spike, length of tubing: approx. 150 cm	0406 3759
Intrafix® Air-matic® RV P, with back-check valve for parallel infusions, silicone-coated spike, length of tubing: approx. 180 cm	0406 3783



A system for safety



Create your own custom workstation  
with fm intensive/fm anaesthesia.  
Just ask – we'd be delighted to advise you!

**B | BRAUN**  
HOSPITAL CARE

B. Braun Melsungen AG  
P.O. Box 11 20  
D-34209 Melsungen  
Phone +49 (56 61) 71-0  
Fax +49 (56 61) 71-3435  
www.bbraun.com