SLE2000 Infant Ventilator



SLE2000 - The Reliable Solution for Infant Ventilation

For many years, the SLE2000 has set the standard in neonatal ventilation.

Using a unique valveless system, the principle of operation of the SLE2000 helps to eliminate inadvertent PEEP, and aid in the total clearance of expired gases.

- ✓ Ability to maintain pressure waveforms at all rates
- ✓ Sensitive airway pressure trigger
- ✓ Comprehensive alarm system
- ✓ Built in oxygen analyser with continuous digital readout
- ✓ Clear and easily set controls
- Excellent pressure square wave means rapid inflation
- √ High pressure dump on alarm
- ✓ Exhaled gases can be filtered
- ✓ Tamper-proof pneumatic controls
- ✓ Optional auxiliary blended flow outlet
- Choice of square or slow rise time (switch selectable)
- Can be used with Nitric Oxide therapy
- √ Easy to clean
- ✓ Low maintenance costs

The smallest thing does matter...

The SLE2000 is specifically designed for neonatal and infant patients. It allows the user to deliver conventional ventilation i.e. CPAP, CMV, PTV or SIMV.

Exceptional pneumatic performance is further enhanced by a patented valveless system.

The valveless system gives the SLE2000 many advantages over conventional ventilators including a less complicated patient circuit with no diaphragms or other obstructions in the exhalation limb of the circuit, helping to reduce expiratory resistance.

All of which makes it easier to assemble and clean and ultimately reduces the potential for errors or wrong assembly.

56 The valveless system makes it easy to clean and assemble 99

Customer Support

Throughout the 50 plus year history of SLE, the company has always put customers and patients first.

We believe that developing life-support equipment should not be taken lightly, and that only the highest quality engineering is acceptable for our products.

Downtime is not desirable in the modern NICU. Logical product design means that routine servicing can be carried out quickly and effectively.

Our teams of qualified engineers support our distributors and end-users world-wide to ensure that your SLE product is always available for optimum performance.

Clinical Education

We know that a ventilator is a vital part of your equipment and that you need to feel comfortable using it. That's why we employ experienced Clinical Specialists in our Clinical Education Department.

These Specialists support our world-wide network of distributors and can help you become more familiar with your SLE products.

Pressure Display

High intensity LEDs show pressure changes.

Mode Selector ···

Switches between modes.

Alarm panel ····

Control and manage all alarms.

Patient connection panel

All patient connections are on the front panel for easy access.

System Pressure controls ...

Large, lockable knobs for easy adjustment.

Oxygen Control ·····

Easy-to-use Oxygen control.



Ordering details

Z2002/FL SLE 2000 with singleuse patient circuits

(International)

Z0002 RS-232 serial port **Z0003** 0-15 l/min blended

output

Z0004/01 Factory fitted

modification to allow BPM range of 126 - 250

BPM

Z0005 External moisture trap

(on air supply only)

N2391 Patient Circuit - Re-usable

(10 mm) (Other types available on request)

N2587/000/050 Bacteria Filter - Single use

(Box of 50)

N2029 Bacteria Filter - Autoclavable N2186 Silencer - Autoclavable

N2190 Exhalation block - Autoclavable

N2002/FL Service manual
N2000/FL User manual
N2001/FL Quick Setup Guide

Compressor

A free-standing air compressor is also available for use with the SLE2000. Please see your local distributor for further details.

BC2188/400/15

Single use breathing circuit for use with SLE2000 and SLE2000 HFO infant ventilators. Temperature port 400 mm from ET manifold (single use). Circuit comes complete with water trap, filter connection kit and adaptors. (Box of 15)

BC6216

Nitric Oxide delivery kit, set of connectors (Paediatric delivery).

BC2508

Nebuliser kit (Paediatric delivery).

BC2110/KIT

Nitric Oxide adaptor kit for BC2188/400 breathing circuits (SLE2000 and SLE2000 HFO infant ventilators). (Box of 5)

N4110

Nitric Oxide dual hose scavenging filter assembly for SLE2000 infant ventilators. (Box of 1)

Ventilator Controls

Modes: Off/Alarm Test, CPAP, CMV, PTV,

SIMV

Breathing Rates: 1 - 125 (126 - 250

with Option Z0004)
Max Inspiratory Time: 0.1 - 3.0 seconds

(0.01 - 0.3 seconds with Option Z0004)

I.E. Ratio: 9.9:1 to 1:99 calculated from

BPM and Ti settings

Pressure Wave: Slow or fast rise,

switch selectable

Manual Breath: Gives single breath in CPAP, CMV, PTV,

SIMV modes Adjustable threshold

PTV Sensitivity: Adjustable threshold Oxygen Blender: 21 - 100% $O_2 \pm 3\%$ Pressure Display

Switch:

CPAP Regulator: Inspiratory Pressure

Regulator: Alarm Mute: Alarm Reset: Max - Mean - Min Range: 0 - 15 mbar 0 - 60 mbar 60 seconds Resets all alarms

(except system fail)

Indicator LEDs

Power: Green LED,

System Fail:

Power on
Red LED,
Main processor
system fail
Trigger Back up:

Green LED,

Indicates a machine delivered breath

Proximal airway LED

Pressure meter: Range -4 to +62

mbar

Digital

BPM Display: Breaths per minute Insp Time Display: Inspiratory time Inspiration/

Expiration ratio Stronger Description (Control of Control of Contr

21 - 100% 0 - 62 mbar

Pressure:

Gauge

Driving Pressure

Gauge: 0 - 60 mbar plus

PEEP level

Alarms

CPAP: Visual and

intermittent audible

PIP / Cycle Fail / Low: Visual and

intermittent audible

Fresh Gas Block or Leak: Visual and

continuous audible

Loss of Mains supply: Battery powered,

audible .

Loss of Air or O₂ supply: Pneumatic, audible

from blender

System Fail: Visual and

intermittent audible

Supplies, Dimensions and Weights

Air and O_2 : 3 - 5 bar

Voltage: 100-120 V 50/60 Hz 220-240 V 50/60 Hz

Power: 120 VA

Fuses: 100-120 V = T500 mA

220-240 V = T200 mA
Protection: Class I Type B

Size, ventilator only: 37 cm W x 31 cm H

x 32 cm D

Height on pole: 137 cm Weight, ventilator only: 10 kg

Designed and manufactured to conform to all relevant international standards for medical devices.

Data subject to change without notice.

SLE is a world leader in the design and manufacture of neonatal ventilation systems.

Years of ventilation experience have given the company an understanding of the challenges facing nurses and clinicians when caring for the tiniest and most critical babies.

From being the pioneers of neonatal Patient Triggered Ventilation (PTV) in the early 1980's, to the introduction of combined HFO (High Frequency Oscillation) in the 1990's, SLE has maintained a position of strength in neonatal ventilation.

The company's guiding principle is to support clinical and nursing staff in their everyday work.

The knowledge and experience gained during years of development is evident in the SLE2000: the result of SLE's ongoing commitments to innovation, competency and care.

