

T-2100 Treadmill

Quiet Efficient yet powerful and durable Stress Test Treadmill

Features

- Long-Surface - Patient can walk on the treadmill without altering their walking style
- Zero Start - A powerful motor is able to increase speed slowly and smoothly
- Both adult and pediatric rails provide support for the user should they require it during use
- Smooth & Quiet - Reduced obtrusive noise
- Self-Calibration - The the treadmill settings never need to be calibrated
- 204 kg (450-pound) weight capacity and speed range of 0-21,7 km/h (0-13.5 miles per hour)

