



DreamStar Auto-CPAP (Non-Evolve)

- **Data** Respiratory events (Apnea, Central Apnea, Hypopnea, Central Hypopnea, Inspiratory Flow Limitation, Snoring), patient flow rate, summary & trends.
- **Data Transfer** : USB Cable, Memory Card
- **Telemonitoring** : Data Transmission via Cellular Modem
- **Humidifier** : Optional (in Evolve Version)
- **Ramp** : Intelligent (I), Time (T)
- **Comfort Mode** : ComfortCalibration Plus (CC+)
- **Tubing** : 22mm, 15mm, Specific
- **Altitude Adjustment** : Automatic
- **Weight** : 1.5 kg (With Power Supply & Humidifier)
- **Country of Origin** : France




DreamStar Auto-CPAP (Evolve)

- **Data:** Respiratory events (Apnea, Central Apnea, Hypopnea, Central Hypopnea, Inspiratory Flow Limitation, Snoring), patient flow rate, summary & trends.
- **Data Transfer :** USB Cable, Memory Card
- **Telemonitoring :** Data Transmission via Cellular Modem
- **Humidifier :** Optional (in Evolve Version)
- **Ramp :** Intelligent (I), Time (T)
- **Comfort Mode :** ComfortCalibration Plus (CC+)
- **Tubing :** 22mm, 15mm, Specific
- **Altitude Adjustment :** Automatic
- **Weight :** 1.5 kg (With Power Supply & Humidifier)
- **Country of Origin :** France



DreamStar Duo (S)

- **Adapt Screen** Helps with adjustment of ventilation parameters (Vt, RR, Avg. Leak, Ti)
- **Modes** : CPAP, & S
- **IPAP** : 3-25 cmH₂O
- **EPAP** : 3-20 cmH₂O
- **Rise Time** : 1-5
- **iSens / ESens** : 1-10
- **Data** : Respiratory Events (Apnea, Hypopnea), Patient Flow, Summary & Trends
- **Data Transfer** : Data Card, USB Cable
- **Telemonitoring** : Data Transmission via Cellular Modem (Optional)
- **Humidifier** : Optional (in Evolve Version)
- **Tubing** : 22mm, 15mm, Specific
- **Altitude Adjustment** : Automatic
- **Weight** : 1.5 kg (With Power Supply & Humidifier)
- **Country of Origin** : France



Dreamstar Duo ST

- **Adapt Screen** : Helps with adjustment of ventilation parameters (Vt, RR, Avg. Leak, Ti)
- **Modes** : CPAP, S & ST
- **IPAP** : 3-25 cmH₂O
- **EPAP** : 3-20 cmH₂O
- **Rise Time** : 1-5
- **iSens / ESens** : 1-10
- **Data** : Respiratory Events (Apnea, Hypopnea), Patient Flow, Summary & Trends
- **Data Transfer** : Data Card, USB Cable
- **Telemonitoring** : Data Transmission via Cellular Modem (Optional)
- **Humidifier** : Optional (in Evolve Version)
- **Tubing** : 22mm, 15mm, Specific
- **Altitude Adjustment** : Automatic
- **Weight** : 1.5 kg (With Power Supply & Humidifier)
- **Country of Origin** : France



EcoStar Auto-CPAP

- **Data** : Respiratory events (Apnea, Central Apnea, Hypopnea, Central Hypopnea, Inspiratory Flow Limitation, Snoring), patient flow rate, summary & trends.
- **Data Transfer** : USB Cable
- **Telemonitoring** : Data Transmission via Cellular Modem
- **Humidifier** : Optional – Compatible with GKH2O
- **Ramp** : Intelligent (I), Time (T)
- **Comfort Mode** : ComfortCalibration Plus (CC+)
- **Tubing** : 22mm, 15mm, Specific
- **Altitude Adjustment** : Automatic
- **Weight** : 750g (Without Power Supply)
- **Country of Origin** : France



EcoStar Info

- **Data** : Pressure usage
- **Humidifier** : Optional
- **Ramp** : Time (T)
- **Patient Circuit** : 22mm
- **Altitude Adjustment** : Auto
- **Country of Origin** : France



EcoStar Intro

- **Humidifier** : Optional
- **Ramp** : Time (T)
- **Patient Circuit** : 22mm
- **Altitude Adjustment** : Auto
- **Country of Origin** : France



S.BOX Auto-CPAP

- **Data** : Respiratory events (Apnea, Central Apnea, Hypopnea, Central Hypopnea, Inspiratory Flow Limitation, Snoring), patient flow rate, summary & trends
- **Data Transfer** : USB Cable, SD-Card, Bluetooth Telemonitoring data transmission via Bluetooth, 3G Modem, Wi-Fi
- **Humidifier** : Optional
- **Ramp** : Intelligent (I), Time (T)
- **Comfort Mode** : Comfort Control Plus (CC+)
- **Tubing** : 22mm, 15mm, Specific
- **Heated Tubing** : 15 mm (Optional)
- **Altitude Adjustment** : Automatic
- **Weight** : 1.3 kg without Power Supply & Humidifier
- **Country of Origin** : France