

CardioMem CM 4000

Simplicity and reliability - two of the outstanding qualities of the CardioMem CM 4000 ECG recorder. This handy, lightweight recorder supports the user throughout all phases of recording and analysing Holter ECGs, from applying the electrodes and entering patient demographic data up to inspecting the ECG waveforms.

Features

- Recording: 2- or 3-channel ECG
- Operating time: Up to 5 days (depending on battery type)
- Sampling rate: 1024 Hz /12 bit
- Display: 2.4" touch screen, self-orientating
- Storage medium: Integrated micro SD card
- Data transfer: Hi-speed USB,
- Casing: Ingress protection against dust & splashed water

Optional

- Voice recorder*
- Bluetooth communication*
- Pacemaker spike detection*
- Thorax impedance channel*
- Activity channel*
- 1024 sample rate

*included as standard on the Advanced recorder.

