

Resting ECG MAC 800

Connecting hearts and minds

MAC 800 is a light, portable ECG with integrated measurements suitable for adults and pediatrics. This 6-channel resting ECG system has the capability to work in both, resting and Arrhythmia mode.

Features

- On-screen 12 SL measurement & interpretation tools
- Bar code and magnetic card reader for patient demographics
- Hook up advisor for optimal ECG signals
- SMS style keypad (alphanumeric) entry
- 6-Channel & full report print-out
- Display: 7" color LCD, 12-lead
- Built-in modem, 2 direct USB for external keyboard and A4 laser print
- Multiple QT correction formulas (Fridericia, Bazett and Framingham)
- R-R analysis to detect underlying HRV patterns and measure R-R intervals
- Long-lasting Li-Ion Battery with 3 hours autonomy

Optional

- Internal (100 or 300 ECGs) & external memory on SD-card
- ECGs can be saved and exported in PDF, XML format
- LAN & Serial port communication to MUSE & CardioSoft for ECG transfer

