

SEER 12

The SEER 12 long-term ECG digital recorder lets you stay connected to your patients for up to 72 hours, so you can better evaluate ongoing cardiac conditions and gather the data you need for real diagnostic confidence. It's another way we connect the hearts of patients to the minds of physicians in ways that may lead to better outcomes.

With portable digital data storage and other advanced features, you can be sure you are getting the high-quality data you need, with a device that gives patients better freedom throughout the process.

Features

- 12-channel ECG with 1,024 Hz sample rate
- Easy-to-read graphic display and user-friendly menu design
- Pacemaker detection
- Saves continuous, high-fidelity ECG data for up to 72 hours
- Displays all 8-leads on-screen for mindful and quality hookups
- Convenient CompactFlash* data storage

