

SunTech Tango+ Blood Pressure

Hands free Monitoring for Stress and Exercise Testing;
designed for Seamless Integration with the GE CASE*
System & CardioSoft* Stress Testing

Features

- Mindful & Reliable: Tolerant to exercise Using SunTech's auscultatory K-sound analysis algorithm
- Seamless Integration: Initiate BP measurements from your stress system & automatically record the data on the printed patient report
- Verify BP Measurements: See & hear the Korotkoff sounds using the on-screen display and the included headphones
- Stat Mode: Continuous BP measurements every 2-10 seconds for time sensitive or emerging situations
- Simple Patient Prep: Simple application with pictorial indicators range and artery indicators
- Choice of Orbit-K* cuff sizes ranging from 18 - 44 cm arm circumference
- Optional SpO₂ module

