

## eBike Ergometers

The connection between GE Healthcare stress-testing systems and eBike allows easy control of the entire examination

### eBike Basic

- Rounds per minute-independent bicycle ergometers
- For load ranges between 20 and 999 watts
- Ergonomic design offers a seating position for optimal performance for all body sizes
- Load accuracy according to DIN VDE 0750-238
- Continuously variable seat height adjustment
- Handlebars adjustment
- Automatic blood pressure measurement (can also be installed later as option)

### eBike Comfort (additional features)

- Electrical, non-graduated seat height adjustment
- Digital display of saddle height
- Dual handlebars adjustment (height and angle)



## eBike EL for Stress Echo Ergometer for load echo cardiography

The eBike EL has a very ergonomic design that provides high patient safety and comfort. The proven eBike intelligence, combined with a very easy to use operator interface, provides clinical focus on the patient during procedures.

### Features

- Programmable and predefined programs
- Electric, variable height adjustment with motor for body height 120-200 cm
- Motorised tilt range from a horizontal position to 45° and lateral from 0° to 45°
- Independent load speed of 20-999 watts
- Large graphic display with load and display of the pulse curve

### Optional

- Blood Pressure Module

